Dr Jean-François Colmas

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Specialist in aesthetic medicine and medical lasers.
Former consultant for the Montpellier medical faculty and
house doctor for the DES at Nice medical faculty.
Graduate of the CIME. Tutor at the International College of
Aesthetic Medicine in Paris. Medical expert and trainer for
Galderma and Fillmed. Laser expert for LUMENIS.
Director of the centre for aesthetic medicine and lasers:
Le Six Foch in Montpellier..

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THE GOLD-STANDARD MEDICAL LASER FOR PROFESSIONALS: SPLENDOR X BY LUMENIS

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Permanent hair removal has been revolutionised thanks to simultaneous beams from two lasers – Alexandrite (755nm)/Nd:YAG (1064nm) – allowing us to increase the Nd:YAG for use up to phototype VI, as well as being an effective treatment for vascular and pigmentation issues.

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The world of laser hair removal has always been divided into two categories of lasers, which are used on the six patient phototypes: phototypes I to III are treated using Alexandrite technology, and phototypes IV to V are treated using the Nd:YAG laser. Phototype III is tricky, as it can be treated using either laser, depending on the laser practitioner and the time of year. Experienced laser practitioners are faced with a difficult choice in terms of the settings for phototypes IIa and IIb, as well as IIIa and IIIb. Today, a revolutionary technological advance has provided a solution: we no longer have to choose between the two categories and can now use the Alexandrite AND the Nd:YAG at the same time and in the same shot. »

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“THE SPLENDOR X IS A LASER CAPABLE OF SATISFYING 100% OF PATIENTS REGARDLESS OF THEIR SKIN COLOUR WITH A VERY CLEAR REDUCTION IN PAIN FOR PATIENTS!”

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Using a double Alex/Nd:YAG shot simultaneously has allowed us to reduce the number of sessions and the pain felt by more than 80%, and to see the disappearance of any contradictory regrowth.

On average, for phototype IIb, we use the 12/7 setting (Alex/Nd:YAG), 6/4 pulse for an aggregate fluence of 19J/cm² in 3Hz from the second session, which gives an acceptable level of pain as well as allowing us to leave 6 to 8 weeks between the second and third sessions, before any partial regrowth.

In my centre, this mixed beam has meant we no longer need to apply a local anaesthetic to the skin prior to the procedure, which allows us to gain time and avoid the risk of burns inherent to an anaesthetic.

This device can be used only with Alexandrite, only with YAG, or with **Alex/Nd:YAG combined**, allowing us to cover all phototypes. It has the added advantage of optimising space in the laser centre. What is more, this YAG laser – which offers sufficient power and is rarely found in the world of lasers – can be used on vascular and pigmentation issues in a quick and effective way.

I believe that the ability to fire 3Hz (3 shots/s), a luxury nowadays, is the gold standard for an experienced hand and enables us to epilate the whole body in just 25 minutes.

We can add two things to this: on the one hand, the revolution of a **square-shaped beam** that offers enough precision for the beard contour as well as for avoiding tattoos or moles, and reduces the risk of overlapping by 90%, something that is inevitable with a round beam; and on the other hand, its comfortable size, which goes up to 27 x 27mm.

In addition to this, it offers important safety features such as reducing the risk of burns and the clear reduction in pain for all patients when using the mixed setting for phototypes I to III. Our centre also noticed a reduction in the

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
number of sessions required, with a very clear efficacy/ satisfaction ratio from the third session, which becomes unanimous from the sixth session.

Focus: A revolutionary asset that is unique on the laser market, the laser's handpiece includes suction with an effective HEPA filter, positioned close to the skin, which eliminates any residual odour and, above all, prevents the laser practitioner from inhaling the smoke released during the shots, which can be damaging for their health. This remarkable feature should be included in the settings of any laser devices designed for hair removal. Another non-negligible financial benefit is that there are no consumables required during the sessions. The flash lamp only needs changing every 1.2 million shots and can be done in-house, thanks to the quick-changing technology.

In conclusion, the future of lasers is undeniably the multiplatform laser, which allows us to treat all phototypes in an optimal way by using a mixed shot for phototypes I to III – a real revolution – while allowing us to treat up to and including phototype VI with the Nd:YAG laser alone. The power of the YAG laser allows us to offer high-quality vascular treatments and treat pigmentation marks with a power reserve that is unparalleled on the market to date. ■